

SÍMBOLOS			
ORI:GEN TECHNOLOGY	HIDROGEL DE SILICONA	HIDROGEL	DISEÑO RX
MENSUAL	TRIMESTRAL	CONVENCIONAL	FILTRO UV
ASFÉRICA/ESFÉRICA	TÓRICA	MULTIFOCAL	MULTIFOCAL TÓRICA

Saphir RX



Mensual
1 blister 3 pack

		MATERIAL	HIDROGEL DE SILICONA, FILCON V 3, 75% H ₂ O (NO IÓNICO)				
		Dk	60				
	DIÁMETRO (mm)	RADIO BASE (mm)	ESFERA (D)	CILINDRO (D)	EJE (°)	ADICIÓN (D)	© DISEÑO
RANGO DE PARÁMETROS	13.00 >> 16.00	6.80 >> 9.80	±30.00	-0.75 >> -8.00	Todos	+0.50 >> +4.00	CD — CN
PASOS	0.50	0.30	0.25	0.25	1º	0.50	—

Gentle80



		MATERIAL	ORI:GEN TECHNOLOGY, 80% H ₂ O				
		Dk	60				
		TINTE DE MANIPULACIÓN	SÍ				
		FILTRO UV	SÍ				
	DIÁMETRO (mm)	RADIO BASE (mm)	ESFERA (D)	CILINDRO (D)	EJE (°)	ADICIÓN (D)	© DISEÑO
RANGO DE PARÁMETROS	13.00 >> 16.00	7.10 >> 9.80	±30.00	-0.75 >> -8.00	Todos	+0.50 >> +4.00	CD — CN
PASOS	0.50	0.30	0.25	0.25	1º	0.50	—

Gentle 59



		MATERIAL	ORI:GEN TECHNOLOGY, 59% H ₂ O				
		Dk	30				
		TINTE DE MANIPULACIÓN	SÍ				
		FILTRO UV	SÍ				
	DIÁMETRO (mm)	RADIO BASE (mm)	ESFERA (D)	CILINDRO (D)	EJE (°)	ADICIÓN (D)	© DISEÑO
RANGO DE PARÁMETROS	13.00 >> 16.00	7.10 >> 9.80	±30.00	-0.75 >> -8.00	Todos	+0.50 >> +4.00	CD — CN
PASOS	0.50	0.30	0.25	0.25	1º	0.50	—

BLU:GEN



		MATERIAL	HIDROGEL DE SILICONA, FILCON V 3, 75% H ₂ O (NO IÓNICO)				
		Dk	60				
		TINTE DE MANIPULACIÓN	SÍ				
		FILTRO UV	SÍ FILTRO SELECTIVO DE LUZ AZUL				
	DIÁMETRO (mm)	RADIO BASE (mm)	ESFERA (D)	CILINDRO (D)	EJE (°)	ADICIÓN (D)	© DISEÑO
RANGO DE PARÁMETROS	11.50 >> 16.50	6.50 >> 9.80	±30.00	-0.75 >> -8.00	Todos	+0.50 >> +4.00	CD — CN
PASOS	0.50	0.30	0.25	0.25	1º	0.50	—



MATERIAL	HIDROGEL DE SILICONA, FILCON V 3, 75% H ₂ O (NO IÓNICO)
----------	--

Dk	60
----	----

	DIÁMETRO (mm)	RADIO BASE (mm)	ESFERA (D)	CILINDRO (D)	EJE (°)	ADICIÓN (D)	☉ DISEÑO
RANGO DE PARÁMETROS	13.00 >> 16.00	6.80 >> 9.80	 ±30.00  ±23.00	-0.75 >> -8.00	Todos	+0.50 >> +4.00	CD — CN
PASOS	0.50	0.30	0.25	0.25	5°	0.50	—

Quattro



MATERIAL	HIDROGEL, FILCON I 1, 49% H ₂ O
----------	--

Dk	15
----	----

TIPO DE MANIPULACIÓN	SÍ
--------------------------------------	--------------------

FILTRO UV	SÍ
-----------	----

	DIÁMETRO (mm)	RADIO BASE* (mm)	ESFERA* (D)	CILINDRO* (D)	EJE (º)	ADICIÓN (D)
RANGO DE PARÁMETROS	 13.00 & 14.50  14.00	 7.70 >>9.80 Ø 13.00 7.10 >>9.20  8.00 >>9.00	 ±30.00  +8.00 >>-12.00	 -0.75 >>-8.00	Todos	A – B – C
PASOS	—	0.30 ( 0.20)**	0.25	0.25	5º	—

Xtensa



MATERIAL	HIDROGEL, FILCON IV 1, 55% H ₂ O
----------	---

Dk 19

TIPO DE MANIPULACIÓN	SÍ
--------------------------------------	--------------------

	DIÁMETRO (mm)	RADIO BASE (mm)	ESFERA (D)	CILINDRO (D)	EJE (°)	ADICIÓN (D)	☉ DISEÑO
RANGO DE PARÁMETROS	☉ ₁ ☉ ₂ ☉ ₃ ☉ ₄ 14.40 ☉ 14.30	☉ ₁ ☉ ₂ ☉ ₃ ☉ ₄ 8.70 ☉ 8.60	☉ +8.00 >>-10.00 ☉ +6.00 >>-8.00 ☉ ₁ ☉ ₂ ±10.00 ☉ ±10.00	☉ -0.75 >>-2.75 ☉ ₁ ☉ ₂ -0.75 >>-5.75	Todos	LOW — HIGH	CD — CN
PASOS	—	—	0.25 & 0.50*	0.50	☉ 10° ☉ ₁ ☉ ₂ 10° & 5°*	—	—

Moldeada 6 pack

Quattro



MATERIAL	HIDROGEL, FILCON I 1, 49% H ₂ O
----------	--

Dk	15
----	----

TIPO DE MANIPULACIÓN	SÍ
--------------------------------------	--------------------

FILTRO UV	SÍ
-----------	----

	DIÁMETRO (mm)	RADIO BASE* (mm)	ESFERA* (D)	CILINDRO* (D)	EJE (º)	ADICIÓN (D)
RANGO DE PARÁMETROS	  13.00 & 14.50  14.00	 7.70 >>9.80 Ø 13.00 7.10 >>9.20  8.00 >>9.00	 ±30.00  +8.00 >>-12.00	 -0.75 >>-8.00	Todos	A – B – C
PASOS	—	0.30 ( 0.20)**	0.25	0.25	5º	—

Anual
1 botella